

NA-Off Reagent

For research use only

Components

Cat. No.	P9121	P9122
NA-Off Reagent	250 ml	500 ml

Storage

Store at room temperature.

Features:

Safe, non - toxic, non - oxidation, non - corrosion, direct contact with the skin; Rapid and effective removal of nucleic acid contamination.

Usage:

- 1) Clean all kinds of table and instrument surface: Spray the product directly on the table and instrument surface, wipe it with ordinary absorbent paper, and dry it.
- 2) Cleaning of pipette: Disassemble the pipette according to the pipette instruction, soak the part of nucleic acid that is liable to contaminate in the product for 5-10 minutes, then thoroughly soak and rinse with distilled water, dry it, and put it back into the pipette.
- 3) Cleaning of glass and plastic containers: Soak the containers in this product, take them out after 5-10minutes, thoroughly soak and rinse them with distilled water, and then dry them for later use.
- 4) Aerosol cleaning: When there is mild aerosol pollution in the laboratory, this product can be used to spray the laboratory comprehensively. After 1-2 hours, wipe the floor and the table with 1:50 diluted 84 disinfectant to effectively remove aerosol pollution; When the laboratory aerosol pollution is serious, 2-3 times of treatment is recommended.