

Super TaqPlus Green PCR Mix

For research use only

Cat. No./Spec.

Cat. No	K034-A	K034-B	K034-C
50- μ l reaction Nos.	40 rxns	200 rxns	4,000 rxns

Description

Super TaqPlus Green PCR Mix is formulated with Invitrogen PlatinumII Taq Hot-Start DNA Polymerase and Platinum SuperFill DNA Polymerase in a certain ratio, and has been verified by NGS to have similar fidelity performance to Platinum SuperFill Green PCR Master Mix. And the amplification performance is stronger. Universal annealing temperatures can be used, which reduces the number of reaction optimization steps and enables simultaneous amplification of different PCR reactions. It has a high capacity for continuous amplification and is more resistant to PCR inhibitors. It also allows for fast-cycling protocols and amplification of long targets (up to 20 kb).

Components

Component	K034-A	K034-B	K034-C
2X Super TaqPlus Green PCR Mix	1 ml	1 ml \times 5	100 ml

This product contains two electrophoresis indicators, blue and yellow, PCR amplification products can be directly electrophoretic.

Features

Superior fidelity, > 150X Taq

60°C universal primer annealing temperature

Excellent specificity, sensitivity, and yield

Reliable amplification of difficult-to-amplify targets (e.g., targets with poor purity, targets with 65% GC content, long PCR requirements)

Storage

Store at -20°C for 2 years.

Protocol

1. Preparation of reaction solution

Add the following reagents to the proper thermal cycler reaction tube or plate on ice or at RT:

Component	50- μ l rxn	Final Conc.
2X Super TaqPlus Green PCR Mix	25 μ l	1X
upstream primer (10 μ M) ^[1]	2 μ l	0.4 μ M

downstream primer (10 μ M) ^[1]	2 μ l	0.4 μ M
template DNA ^[2]	1-4 μ l	<0.5 μ g
Water, nuclease-free	to 50 μ l	–

[1] Recommended range of final primer concentration: 0.1-1 μ M. The concentration can be reduced when the specificity is poor, and the concentration can be increased when the efficiency is low.

[2] The optimal dosage varies with different templates. The recommended dosage for some DNA templates is as follows (50 μ l reaction system).

Template	Human genomic DNA	λ DNA	cDNA	Plasmid DNA
Dosage	1ng-500g	0.5ng-5ng	1-5 μ l	0.1ng-10ng

2. Perform PCR using the following thermal cycling condition

Stage	Temperature	Time	Number of Cycles
Initial Denaturation	94°C	2 min	1
Denaturation	94°C	30 sec	25-35
Annealing	60°C ^[1]	30 sec	
Extension	72°C	20 sec/kb	
Final Extension	72°C	5 min	1

[1] 60° C can meet the annealing of most primers, and the appropriate annealing temperature can be found by gradient PCR for special primers.

Quality Control

The absence of endodeoxyribonucleases, exodeoxyribonucleases and ribonucleases is confirmed by appropriate quality tests. Functionally tested in amplification of a single-copy gene from human genomic DNA.

Product Use Limitations

This product is sold exclusively for research purposes and in vitro use. Neither the product, nor any individual components, was tested for use in diagnostic applications or for drug development, nor is it suitable for administration to humans or animals. Please refer to the MSDS, available upon request.