

## RNase H

### Instruction for Use

Cat. No./Spec.: E1016-A/100 U, E1016-B/500 U

Concentration: 5 U/ $\mu$ L

#### Product Description

Ribonuclease H (RNase H) can specifically degrade the RNA strand within RNA-DNA hybrids. This enzyme does not hydrolyze the phosphodiester bonds in single-stranded or double-stranded DNA and RNA.

#### Components

Component	E1016-A	E1016-B
RNase H (5 U/ $\mu$ L)	20 $\mu$ L	100 $\mu$ L
10X RNase H Buffer	1 mL	1 mL

#### Storage Condition & Shelf Life

Store at -20°C.

#### Source

*E. coli* strain MRE-600.

#### Unit Definition

A unit is defined as the amount of enzyme required to catalyze the formation of 1 nmol of acid-soluble products within 20 minutes at 37° C.

Enzyme activity is determined in the following mixture: 20 mM Tris-HCl (pH 7.8), 40 mM KCl, 8 mM MgCl<sub>2</sub>, 1 mM DTT, 24  $\mu$ M [<sup>3</sup>H]-poly(A) ·poly(dT), 0.03 mg/mL BSA, and 4% (v/v) glycerol.

#### Scope of Application

- Removal of mRNA before the synthesis of the second-strand cDNA

- RT-PCR and qRT-PCR: Removal of RNA after the first-strand cDNA synthesis
- Removal of poly(A) sequence from mRNA after hybridization with Oligo(dT)
- Site-specific cleavage of RNA
- Study of products from *in vitro* polyadenylation reactions

#### Inhibition and Inactivation

- Inhibitors: metal chelating agents, SH-blocking reagents.
- Heat at 65°C for 10 minutes.

#### Protocol

Second-strand cDNA synthesis

1. Conduct the first-strand cDNA synthesis reaction according to the recommendations for the specific reverse transcriptase being used.
2. Add the following components (on ice) to the 20  $\mu$ L first-strand cDNA synthesis reaction mixture:

Component	Amount
10X reaction buffer for DNA Polymerase I*	8 $\mu$ L
RNase H (5 U/ $\mu$ L)	0.2 $\mu$ L
DNA Polymerase I	30 U
ddH <sub>2</sub> O	To 100 $\mu$ L

\* 10X reaction buffer for DNA Polymerase I: 500 mM Tris-HCl (pH 7.5 at 25 °C), 100 mM MgCl<sub>2</sub>, 10 mM DTT.

- ③ Gently mix and briefly centrifuge.
- ④ Incubate at 15° C for 1 hour. Do not exceed the temperature of 15° C.
- ⑤ Add 2.5  $\mu$ L (12.5 U) of T4 DNA Polymerase (#K011) and incubate at 15° C for 5 minutes.
- ⑥ Terminate the reaction by adding 5  $\mu$ L of 0.5 M EDTA (pH 8.0). The blunt-end cDNA purified by phenol/chloroform can be used for further cloning-related procedures, such as adapter ligation, phosphorylation, size separation, ligation, and transformation.

This product is for research use only.