

DNase I, RNase-Free, HC

Instruction for Use

Cat. No./Spec.: E1018-A/1,000 U

Concentration: 50 U/ μ L

Product Description

DNase I (without RNase) is an endonuclease that can cleave both single-stranded and double-stranded DNA. It hydrolyzes phosphodiester bonds, producing single deoxyribonucleotides and oligodeoxyribonucleotides with 5'-phosphate groups and 3'-OH groups.

The activity of this enzyme is strictly dependent on Ca^{2+} and is activated by Mg^{2+} or Mn^{2+} ions. In the presence of Mg^{2+} , DNase I cuts each strand of dsDNA in a statistically random and independent manner. When Mn^{2+} is present, the enzyme almost cuts both DNA strands at the same site, resulting in DNA fragments with blunt ends or overhangs of one or a few nucleotides.

Components

Component	E1018-A
DNase I, RNase-Free, HC (50 U/ μ L)	20 μ L
10X DNase I Buffer (Mg^{2+} Plus)	1 mL
EDTA (50 mM)	1 mL

Storage Condition & Shelf Life

Store at $-20^{\circ}C$.

Source

Recombinant *E. coli* strain containing the cloned gene encoding bovine DNase I.

Unit Definition

A unit is defined as the amount of enzyme required to completely degrade 1 μ g of plasmid DNA within 1 minute at $37^{\circ} C$.

Enzyme activity is determined in the following mixture: 10 mM Tris-HCl (pH 7.5 at $25^{\circ} C$), 2.5 mM $MgCl_2$, 0.1 mM $CaCl_2$, and 1 μ g of pUC19 DNA.

One unit of DNase I is equivalent to 0.3 Kunitz units.

Features

- Recombinant enzyme
- Purified from a non-animal host, with low levels of endogenous RNase

Scope of Application

- For the preparation of DNA-free RNA.
- To remove template DNA after *in vitro* transcription.
- To prepare DNA-free RNA before RT-PCR and RT-qPCR.
- In conjunction with DNA Polymerase I for DNA labeling through nick translation.
- For the study of DNA-protein interactions using DNase I (RNase-free) footprinting.
- To generate a library of randomly sheared DNA insert fragments. Reaction buffer containing Mn^{2+} is used.

Inhibition and Inactivation

- *Inhibitors: metal chelating agents, millimolar concentrations of transition metals (such as zinc), SDS (even at concentrations below 0.1%), reducing agents (DTT and β -mercaptoethanol), and ionic strength higher than 50-100 mM.*
- *Inactivation is achieved in the presence of EGTA or EDTA by heating at $65^{\circ}C$ for 10 minutes (using at least 1 mol of EGTA/EDTA per 1 mol of Mn^{2+}/Mg^{2+}).*

Protocol

1. Removal of genomic DNA from RNA samples

① Prepare the reaction system in an RNase-Free tube:

Component	Amount
RNA	1 μ g

10X DNase I Buffer (Mg ²⁺ Plus)	1 μ L
DNase I, RNase-Free, HC (50 U/ μ L)	1 μ L
DEPC-treated Water (#R2041)	To 10 μ L

② Incubate at 37°C for 30 minutes.

③ Add 1 μ L of 50 mM EDTA and incubate at 65°C for 10 minutes. In the absence of chelating agents, RNA is hydrolyzed with divalent cations during heating. Alternatively, use phenol/chloroform extraction.

④ Use the prepared RNA as a template for reverse transcriptase.

Note:

- Use no more than 1 U of DNase I, RNase-Free, per 1 μ g of RNA.

- DNase I, RNase-Free, HC can be diluted in 1X reaction buffer or in storage buffer before use to extend storage time.

- The volume of the reaction mixture and the 50 mM EDTA solution can be scaled up to obtain larger quantities of RNA. The recommended final concentration of RNA is 0.1 μ g/ μ L.

- The RNase inhibitor RNasin (#R2011) can also be included in the reaction mixture to prevent RNA degradation.

2. Removal of template DNA after in vitro transcription

① Add 2 U of DNase I, RNase-Free, per 1 μ g of template DNA directly to the transcription reaction mixture. In some cases, the amount of enzyme should be determined empirically.

② Incubate at 37°C for 15 minutes.

③ Inactivate DNase I by phenol/chloroform extraction.

3. DNA labeling by nick-translation

① Prepare the reaction system:

Component	Amount
10X reaction buffer for DNA Polymerase I	2.5 μ L
Mixture of 3 dNTPs, 1 mM* each (without the labeled dNTP)	1.25 μ L
[α - ³² P]-dNTP,	1.85-3.7 MBq

~110 TBq/mmol (3000 Ci/mmol)	(50-100 μ Ci)
DNase I, RNase-free	1 μ L
freshly diluted to 0.002 U/ μ L**	
DNA Polymerase I (#E1021)	5-15 U
Template DNA	0.25 μ g
Water, nuclease-free (#P9021)	To 25 μ L

Note:

* To prepare a mixture of three unlabeled dNTPs (each at 1 mM), mix 1 μ L of each dNTP sample (100 mM) with 97 μ L of nuclease-free water. Store at -20°C.

** DNase I, RNase-Free can be diluted with the 1X reaction buffer of DNA Polymerase I: 50 mM Tris-HCl (pH 7.5 at 25°C), 10 mM MgCl₂, and 1 mM DTT.

② Incubate immediately at 15°C for 15-60 minutes.

③ Terminate the reaction by adding 1 μ L of 0.5M EDTA, pH 8.0.

④ Take an aliquot (1 μ L) to measure the efficiency of label incorporation. The expected specific activity of DNA should be at least 108 cpm/ μ g.

This product is for research use only.