

PCR Enhancer

#P9041, 500 μ l

Store at -20°C

For research use only.

In total 1 vial.

Description

PCR Enhancer System is used in conjunction with DNA polymerases to optimize PCR of problematic and/or GC-rich templates. The system offers higher primer specificity, broader magnesium concentration, broader annealing temperature, and improved thermostability DNA polymerases.

Applications

- Complex template with GC-rich
- High efficiency, quality and yield

Usage

Add 2 – 8 μ l PCR Enhancer to 25 μ l amplification system, the other amplification conditions and procedure remain the same.

PRODUCT USE LIMITATION.

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.